HaiGene SYBR Green miRNA定量解决方案

—miRNA提取、反转录、RealTime PCR

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
- Stage 1: 95°C 15min
- Stage 2: 95°C 5s
 55°C 5s
 70°C 30s 30~40 cycles
- Stage 3: Dissociation analysis

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- [1] miR-Q: a novel quantitative RT-PCR approach for the expression profiling of small RNA molecules such as miRNAs in a complex sample. *BMC Mol Biol.* 2008. PMC2374797.
- [2] Specific and sensitive quantitative RT-PCR of miRNAs with DNA primers. *BMC Biotechnol.* 2011. PMC3135530.
- [3] Reliable reference miRNAs for quantitative gene expression analysis of stress responses in *Caenorhabditis elegans*. *BMC Genomics.* 2014. PMID:24656064.
- [4] Profiling microRNAs in lung tissue from pigs infected with *Actinobacillus pleuropneumoniae*. *BMC Genomics.* 2012. PMID:22953717.
- [5] The role of viral and host microRNAs in the Aujeszky's disease virus during the infection process. *PLoS One.* 2014. PMID:24475202.
- [6] MicroRNA polymorphisms and environmental smoke exposure as risk factors for oesophageal squamous cell carcinoma. *PLoS One.* 2013. PMID:24205249.
- [7] Plasma proANP and SDMA and microRNAs are associated with chronic mitral regurgitation in a pig model. *Endocr Connect.* 2013. PMID:24029364.

- [8] Advances in microRNA experimental approaches to study physiological regulation of gene products implicated in CNS disorders. *Exp Neurol*. 2012. PMID:22245616.
- [9] Maternally deposited germline piRNAs silence the tirant retrotransposon in somatic cells. *EMBO Rep*. 2013. PMID:23559065.
- [10] miRNA expression profile analysis in kidney of different porcine breeds. *PLoS One*. 2013. PMID:23372853
- [11] Determination of reference microRNAs for relative quantification in porcine tissues. *PLoS One*. 2012. PMID:22970213.
- [12] Identification of circulating miRNA biomarkers based on global quantitative real-time PCR profiling. *J Anim Sci Biotechnol*. 2012 PMID:22958414.

Material and Methods for miRNA extraction and qPCR

The miRNA was extracted from the tissue or cell samples using miRNA Extraction Kit (HaiGene, Harbin, China) according to manufacturer's protocol. To investigate the miR-21 expression level in tissue samples, cDNA synthesis was carried out according to the protocol of One Step miRNA cDNA Synthesis Kit(HaiGene, Harbin, China), in which poly(A) tailing of the miRNAs is followed by reverse transcription with a tagged poly(T) primer described by Balcells et al [2]. Briefly, the reaction mix (20 μ l) consisting of 100 ng of miRNA, 5 μ l of 4 \times One step miRNA RT Solution and 2 μ l of 10 \times miRNA RT Primer was incubated at 37 $^{\circ}$ C for 1 hour followed by enzyme inactivation at 95 $^{\circ}$ C for 5 minutes. Quantification of miRNA was performed by HG miRNA SYBR Green PCR Kit(HaiGene, Harbin, China) . RealTime PCR was performed in 20 μ l with 1 μ l of cDNA, 1 of each primer, and 4 μ l of Golden HS SYBR Green qPCR Mix under the following conditions:95 $^{\circ}$ C for 15min, 35 cycles of 95 $^{\circ}$ C for 5s, and 60 $^{\circ}$ C for 30s in LightCycler 480 (Roche). MiR-21 levels were normalized to U6B RNA levels using the $2^{-\Delta Ct}$ model. Real-time PCR assays were performed in duplicate for each sample, and the mean value was used for the calculation of miRNA expression levels. The statistical analyses were done using Microsoft Excel (Microsoft) both to calculate the SD and to test for statistically significant differences between the samples using a T test. A value $P < 0.05$ was considered statistically significant.